

GLIDE GEL CONDUCTIVE GEL

NeoGenesis Glide Gel is a safe and highly effective conductive gel that provides easy glide for professional treatments and at-home devices. Pair with NeoGenesis serums for maximum effectiveness.

SIZES: 30mL, 120mL

(size offerings may vary)

PRODUCT BENEFITS

- Ensures safe and effective use of devices
- Enhances delivery of skin enhancing ingredients
- For all skin types
- Allows for customization with the Moisturizing Mist for a thinner viscosity

DIRECTIONS

After cleansing, apply a NeoGenesis serum of choice. Brush the Glide Gel directly onto the skin and follow device directions. Reapply or reactivate with NeoGenesis Moisturizing Mist if gel begins to dry during the treatment. Remove remaining gel with a warm wet cloth when the treatment is complete. Finish with a NeoGenesis moisturizer.

Keep out of eyes. For external use only.

INGREDIENTS

Water, Glycerin, Propanediol, Hydroxyethylcellulose, Benzyl Alcohol, Dehydroacetic Acid, Lactic Acid.

